

## **Finer Face Micro-Current Unit - Frequently Asked Questions.**

Do I have to be a qualified beauty therapist to operate the Finer Face Unit?

No. When you purchase the Finer Face Unit, we fully train you on how to use your new machine to achieve the ultimate results. This gives you the opportunity to start up your own business from home, or in commercial premises.

### **How does the Finer Face Unit work?**

Finer Face uses totally safe and painless low-level electrical impulses to re-programme the muscles of your face, to revitalise these muscles and to assist your skin tissues in becoming healthier as your blood circulation is improved. Your body is like a machine run by electrical signals from your heart and your thoughts. Even the tone of your skin and muscles are related to the body's electrical system. The skin's ability to maintain its tone and elasticity is controlled by chemical reactions within and between the cells of your skin, and these reactions are triggered by electrical impulses.

As the body ages these reactions slow down and the skin loses its elasticity and suppleness. The underlying muscles, which are the skin's support, also begin to lose their tone, and this can result in the apparent sagging of the skin. The Finer Face Unit uses very low variable electronic currents to stimulate your body's tiny currents which result in a visible firming and smoothing of the skin. The system uses two twin-tipped probes to gently stimulate the muscles in your face and neck. There is no pain or discomfort. You will see the results immediately and your skin will glow and appear to keep lifting for at least 24 hours after the treatment has been performed.

### **How many treatments will my clients require?**

For the best, long-term results it is highly advisable for clients to complete a series of 12 treatments over a period of 6 weeks. This requires following with a maintenance treatment plan every month thereafter. It is vital to complete 12 treatments as they are made up of 4 low, 4 medium, and then 4 high intensity treatments.

### **Does the Finer Face Unit work off the mains power?**

No. The Finer Face Unit is powered by 6 D Cell batteries which provide you with the necessary micro-current required to perform all of the pre-set programmes specifically designed for the face. These 6 batteries will provide you with 250 hours of treatment time. The Finer Face Unit will not permit you to do a treatment if the batteries are too low, as the LCD screen will tell you when replacement is necessary.

## **Why is the Finer Face Unit battery operated, and not mains supplied?**

Micro-current or MENS (Micro-Current Neuromuscular Stimulation) is a similar current to which our body produces. If a current is too strong or if the facial muscles are physically stimulated with small electrodes which are usually attached directly to the skin, the treatment will drop out. MENS is one millionth of an amp, (for example it would take 1.5 million micro-current units to light up a 60-watt bulb). The Finer Face Unit is designed with the correct percentage of micro-current to ensure acceptance from the tiny facial muscles. Your client will receive a relaxing treatment with long lasting results, rather than with a superficial difference that only lasts for a few hours.

## **Does the Finer Face Unit only have 1 wave shape like most of the other MENS units?**

The Finer Face Unit has 4 wave forms which vary through the different treatments. These wave forms are varied to penetrate skin tissue to a pre-determined depth which allows the machine to isolate specific tissues for re-programming. The low current stimulates the body's own tiny currents and also avoids any unpleasant stinging that other machines sometimes inflict.

## **How do I know which wave form to choose?**

The Finer Face Unit eliminates any guess work as it is fully computerised, the correct wave form is pre-programmed into each individual treatment so it is automatically selected for you.

## **How many different treatments are built in to the Finer Face Unit?**

The Finer Face Unit has 9 individual pre-set treatments to choose from. The in-built computer automatically adjusts to the frequency, rate and wave form for whichever treatment you choose.

## **Why are dual tipped probes better than single tips?**

The dual tipped probes on the Finer Face Unit gives you an interferential cross current which single tipped probes cannot do. It will also cut the treatment time down.

## **Does the Finer Face Unit check the resistance of the skin?**

The Finer Face Unit checks the skins resistance 1024 times per second. This ensures automatic correct current depth for whichever treatment you are performing.

**What type of treatments can I do with the Finer Face Unit?**

Non-Surgical Face Lifting, Acne Expulsion, Pigmentation Normalisation, Deep Line Softening, Fine Line Depletion, Exfoliation, Disencrustation, Ionization, Scar Tissue Repair, Couperose, Capillaries, Rosacea, Eczema, Psoriasis, Dermatitis, Cuts, Burns & Stretch Marks, Hand Rejuvenation & Décolleté.

**What makes the Finer Face Unit stand out against its competitors?**

The Finer Face Unit is a multi-unit. It has been designed to ensure that all problematic areas on the face can be addressed. The Finer Face Unit has been developed by a leading Micro-Current Institute and tested on over 2000 people, with a 96% success rate after just one treatment.